



**Demand a Higher Standard**

# Centrifuge Comparison Chart



Tel: (713) 926-2623  
Fax: (713) 926-7736

Petroleum Testing Equipment  
Suppliers Worldwide Since 1930

<u>Product Feature</u>	<u>Melton</u>	<u>Transport Series</u>	<u>Benchmark 2000</u>	<u>Benchmark S</u>	<u>Benchmark C</u>
Tube Style Compatibility	B, C, D, E	A, B, D	B, C, D, E	B, C, D, E	B, C, D, E
Voltage	120/240VAC	12VDC or 115VAC	120/240VAC	120/240VAC	120/240VAC
Max RPM	1750	1800	1800	2000	2300
Manual RPM/Speed Adjustment		✓			
Manual RPM/Speed Adjustment with LCD Display		9100 only	✓	✓	✓
Heated Operation	✓	✓	✓	✓	✓
Tube Preheater Pocket		✓			
High Intensity Heat (Dual Heating Elements)					✓
Vent Chamber				✓	✓
Insulated		✓	✓	✓	✓
Manual Temperature Adjustment with LCD Display (°F only)		9100 only	✓	✓	
Digital Programmable Temperature Controller with LCD Display					✓
Digital Programmable Countdown Timer with LCD Display			✓	✓	✓
Rated for Outdoor Use	✓	✓			
Class I, Division 2 Certification					✓
Emergency Stop Button			✓	✓	✓
Automatic Lid-Locking during Operation			✓	✓	✓
Plug-in (Power Cord) Option			✓	✓	✓