



Demand a Higher Standard



Tel: (713) 926-2623
Fax: (713) 926-7736

*Petroleum Testing Equipment
Suppliers Worldwide Since 1930*

For Safety's Sake: Why Class I, Division 2 Certification Matters

A Class I, Division 2 location is deemed hazardous due to the presence of flammable gases – potentially in quantities high enough to produce an explosive or ignitable mixture. A classified hazard does not normally exist in that location, but can possibly appear under abnormal conditions.

The ETL (Electrical Testing Laboratories) Mark on L-K Industries products indicates certification was completed by Intertek, a third-party NRTL (Nationally Recognized Testing Laboratory). All products certified for use in Class I, Division 2 sites meet the applicable published standards and are safe in these potentially hazardous locations.

Class I, Division 2 Certified Products

	Product Name	Model #	Rating
Centrifuge	Benchmark C Laboratory Centrifuge (C)	LB14-C120C	120 VAC, 50/60 Hz, 505W
	Benchmark C Laboratory Centrifuge (C)	LB14-C240C	240 VAC, 50/60 Hz, 505W
	Benchmark C Laboratory Centrifuge (D)	LB14-D120C	120 VAC, 50/60 Hz, 505W
	Benchmark C Laboratory Centrifuge (D)	LB14-D240C	240 VAC, 50/60 Hz, 505W
Sample Heater	Criterion Water Bath (1 Gal)	WB-T2-110	110 VAC, 50/60 Hz, 900W
	Criterion Water Bath (1 Gal)	WB-T2-220	220 VAC, 50/60 Hz, 900W
	Criterion Water Bath (C)	WB-T14C-110	110 VAC, 50/60 Hz, 900W
	Criterion Water Bath (C)	WB-T14C-220	220 VAC, 50/60 Hz, 900W
	Criterion Water Bath (D)	WB-T12D-110	110 VAC, 50/60 Hz, 900W
	Criterion Water Bath (D)	WB-T12D-220	220 VAC, 50/60 Hz, 900W
	TradeMark Heater (A)	LAB-H-A-110	110 VAC, 50/60 Hz, 400W
	TradeMark Heater (A)	LAB-H-A-220	220 VAC, 50/60 Hz, 400W
	TradeMark Heater (B)	LAB-H-B-110	110 VAC, 50/60 Hz, 400W
	TradeMark Heater (B)	LAB-H-B-220	220 VAC, 50/60 Hz, 400W
	TradeMark Heater (C)	LAB-H-C-110	110 VAC, 50/60 Hz, 400W
	TradeMark Heater (C)	LAB-H-C-220	220 VAC, 50/60 Hz, 400W
	TradeMark Heater (D)	LAB-H-D-110	110 VAC, 50/60 Hz, 400W
	TradeMark Heater (D)	LAB-H-D-220	220 VAC, 50/60 Hz, 400W
	Refined Fuel Test Cylinder Heater	LAB-H-S-110	110 VAC, 50/60 Hz, 400W
	Refined Fuel Test Cylinder Heater	LAB-H-S-220	220 VAC, 50/60 Hz, 400W
	LP Gas Corrosion Test Cylinder Heater	E-4-110-D	110 VAC, 50/60 Hz, 650W
	LP Gas Corrosion Test Cylinder Heater	E-4-220-D	220 VAC, 50/60 Hz, 650W